



DIFFERENTIAL WORKSHEET 2010

The slides have been examined and found to be homogeneous for the purpose of these exercises.

JANUARY	Cycle 10 Run 1	BF10-01	Differential Count
FEBRUARY	Cycle 10 Run 2	BF10-02/MO10-02a	Differential Count
MARCH	Cycle 10 Run 3	BF10-03	Differential Count
APRIL	Cycle 10 Run 4	BF10-04	Differential Count
MAY	Cycle 10 Run 5	BF10-05/MO10-05a	Differential Count
JUNE	Cycle 10 Run 6	BF10-06	Differential Count
JULY	Cycle 10 Run 7	BF10-07/MO10-07a	Differential Count
AUGUST	Cycle 10 Run 8	BF10-08	Differential Count
SEPTEMBER	Cycle 10 Run 9	BF10-09	Differential Count
OCTOBER	Cycle 10 Run 10	BF10-10/MO10-10a	Differential Count
NOVEMBER	Cycle 10 Run 11	BF10-11	Differential Count
DECEMBER	Cycle 10 Run 12	BF10-12	Differential Count

- Please perform a differential count on the stained blood film provided.
- The RCPA Haematology QAP currently use the ICSH stain. The stain will be noted on the insert included with each dispatch.
- For the sake of uniformity, please **do not** report bands but include them as either metamyelocytes or neutrophils according to the stage of maturation. **Please report your differential in whole numbers only.**
- When entering results online you must enter a numerical value in all cell lines, ie. if **no** basophils, myelocytes, blasts etc are counted on the film "0" must be entered in the entry field. If a field is left blank it will appear that no results were returned for that cell line.